

## Simply Seamless DNA Assembly Protocol

1. Set up the assembly reaction in a microtube on ice: (keep all components on ice and transfer stocks back to freezers promptly after use)

|                                  | 2-3 Fragment Assembly  | 4-6 Fragment Assembly  | (+) control     |
|----------------------------------|--|--|-----------------|
| Amount of DNA fragments          | 1:1 molar ratio of fragments, 50-200 ng DNA per reaction, volume 8 $\mu$ l or less | 1:1 molar ratio of fragments, 50-200 ng DNA per reaction, volume 8 $\mu$ l or less | 8 $\mu$ l (S95) |
| Simply Seamless Enzyme Mix (S92) | 1 $\mu$ l  | 1 $\mu$ l  | 1 $\mu$ l       |
| 10X Reaction Buffer (S94)        | 1 $\mu$ l  | 1 $\mu$ l  | 1 $\mu$ l       |
| H2O                              | Bring volume up to 10 $\mu$ l  | Bring volume up to 10 $\mu$ l  | -               |
| Total Volume                     | 10 $\mu$ l   | 10 $\mu$ l   | 10 $\mu$ l      |

- Incubate the samples in a thermocycler, water bath or a hot block at 37°C for 15-30 minutes. (incubation time may be extended to 60 minutes when more than 4 fragments are being assembled)
- Transfer the completed reaction to ice and use immediately or store at -20°C for subsequent transformation.
- Transform Competent *E. coli* cells (DH5 $\alpha$  provided in the cloning kit or purchased separately, S96) with 3  $\mu$ l of the assembled product, following the transformation protocol.